

FLUORESCENT CELLULAR MODEL FOR ALS/FTD DRUG SCREENING TARGETING TDP-25 AGGREGATES

Patricia Villacé, Rosa M. Mella, Meritxell Roura-Ferrer, Clarisa Salado


Abstract

Amyotrophic lateral sclerosis (ALS) and frontotemporal dementia (FTD) share pathogenic mechanisms, including TDP-25 aggregation, which contributes to neuronal toxicity. We developed a fluorescent cellular model to visualize and quantify TDP-25 aggregates in real time, providing a powerful tool for drug screening. This system enables efficient evaluation of compounds targeting TDP-25 aggregation, offering a promising platform for identifying new ALS/FTD therapies.

Introduction

The accumulation of misfolded protein aggregates is a hallmark of several neurodegenerative diseases, including ALS and FTD. In particular, TDP-43 and its truncated fragment TDP-25 have been implicated in neuronal toxicity in these disorders. Due to the lack of effective treatments, it is crucial to develop cellular models that enable the evaluation of potential drugs aimed at reducing these aggregates.

In this study, we present a fluorescence-based cellular model for real-time detection of TDP-25 aggregates. By expressing a fluorescently tagged construct, we can monitor aggregate formation dynamics and assess the impact of different compounds on their reduction. This approach represents a significant advancement in drug screening, providing an efficient and quantifiable platform for identifying novel therapeutic strategies against ALS and FTD.




Materials and Methods

Cell Line Development: U2OS cells were transfected with TDP25-tGFP cDNA.

Aggregation Assay: Cells were treated overnight with 50 μ M Sodium Arsenite in OptiMEM.

HTS: Cells were incubated overnight with 10 μ M Prestwick library compounds diluted in 50 μ M Sodium Arsenite-OptiMEM.

Imaging: After treatments, PBS replaced the medium images were captured and analyzed using the Thermo Fisher Cell Insight CX7 (Fig. 5).



Results

Stress Induction







Figure 4. Immunofluorescence analysis of TDP-25 localization. The first column shows green fluorescent TDP-25 cells, the second column displays antibody labeling (TDP-43 Antibody), and the third column presents the merged images. Control cells (top row) are compared to arsenite-treated cells (bottom row) to assess changes in TDP25 distribution.



High-ThroughPut Assay



Conclusions

Fluorescence-Based Model Validation: Our model effectively tracked TDP-25 aggregate formation and clearance, proving its value for drug screening in neurodegenerative disease research.

Advancement in Drug Discovery: This study supports the development of cellular models for high-throughput screening, advancing research for TDP-43-related disease therapies.

Compound Efficacy: Two compounds from the library reduced TDP-25 aggregation in a dose-dependent manner, highlighting their potential as therapeutic candidates for ALS and FTD.